

Product Data Sheet

Compressed Gauze, Z-fold

Description: White gauze bandage made of 100 % cotton, 11.43 cm x 3.75 m (4.5" x 4.1 yds), 6-ply, compressed and vacuum-packed, sterile, high absorbency, double-layered packaging for easy and hygienic handling, wide range of applications for wound care.

REF: NR70-VSCGF-01



Product features		Packaging properties, storage	
Material	100 % cotton	Basic package	Content: 1 piece Double-layer packaging made of PE film (outside) and nylon film (inside)
Color	White	Sterile product	Yes
Dimensions (basic package)	Approx. 15 x 10 x 1.5 cm	Shelf life	7 years from date of manufacture
Weight (incl. packaging)	Approx. 40 g	Collective pack	/
Regulatory and market-specific information		Shipping packaging	200 pieces
Regulation	CE-compliant medical device class Is according to Regulation (EU) 2017/745 (MDR)	Shipping packaging: Dimensions, weight	Approx. 54 x 34 x 20 cm, approx. 9 kg
NSN	6510-01-605-2234	Storage information	Protect from sunlight Store in a dry place
Country of manufacture	CN		
Customs tariff number	3005901000		



European Importer:
nextRESCUE GmbH
Brennhorster Str. 12
32479 Hille, Germany



Caretechion GmbH
Niederrheinstr. 71
40474 Düsseldorf
Germany



Anji Hengfeng Sanitary Material Co., Ltd.
Ancheng, Dipu Town, Anji County
313300 Zhejiang, P.R. China